

# Robot Control and Plan Execution for Automated Experiments in Industrial Laboratory

Master Thesis in Artificial Intelligence and Robotics

Department of Computer, Control and Management Engineering  
DIAG "Antonio Ruberti"



**SAPIENZA**  
UNIVERSITÀ DI ROMA

Candidate: Franci Rrapi, 1660405

Thesis Advisor: Prof. Luca Iocchi, Sapienza – University of Rome

Co-Advisor: Dr. Fabio Zonfrilli, Procter & Gamble – Belgium

Academic Year: 2021/2022

# Outline



- Introduction of AIPlan4EU project
- Problem summary
- Software and Tools
- Laboratory and Devices
- Robot configuration
- How to avoid collision
- How the robot grasps the pouch
- Overview of robot movements
- Actions and Fluents implemented
- Plan implementation
- Video showing complete experiments

- AIPlan4EU aims to bring the most advanced AI planning technology to companies
- “Use-cases” drive the design of AI planning systems
- For this thesis, use-case consists of the executive procedure for the analysis of the laundry detergent pouch, performed by a robotic arm, in a dedicated P&G laboratory
- Furthermore, to lay the foundations for an AI planning framework that leverages on the “human-in-the-loop” concept
- To accomplish this, all robot components and devices in the laboratory have been connected between them and designed AI plans

# Problem overview




The **executive procedure for the analysis** of the laundry detergent **pouch** consists in:


- Detect the pouch position
- Grasp the pouch
- Execute the weight test
- Execute the tightness test
- Trash the pouch

We want to **REPEAT** this procedure ***for each*** **pouch** in the drawer and ***for every*** available **drawer**.

# Software and Tools

  
Robot Operating System  
A collection of software libraries and tools that aim to build robot applications

  
GAZEBO  
3D robotics simulator that consent to create realistic and dynamic simulations of robots in a work environment



**docker**  
A state-of-the-art containerization technology that guarantees portability, autonomy, accelerates and simplifies application development

 **MoveIt**  
A robotic manipulation platform that incorporates the latest advances in motion planning, manipulation, 3D perception, kinematics, control and navigation

**YOLOv5**  
One of the most famous object detection tools based on CNN due to its speed and accuracy

# Laboratory and Devices



# Robot configuration

## ***URDF - Unified Robotic Description Format***


XML format to describe the complex interconnected structure of the robot (link in connected in Parent-Child relationship)



New Movelt  
configuration





Inverse Kinematics  
solution




 **Movelt**

# Collision avoidance





# Grasping Pouch



**Note!** If the robot fails to grasp the pouch, it enters in **Recovery Pouch Mode** until he grasp it

# Actions & Fluents

- **GOTO** → performs movements of the robot for reaching various joint states or pose goals that are reachable in the environment (home position “*vision2*”, scale, mark10, bin, etc.)
- **BOX** → performs lots of complex movements of the robot for opening or closing a certain drawer
- **MOVEGRIPPER** → this action is useful for controlling and sending commands to the gripper (activation, reset, open or close gripper, set the speed or force movement)
- **PERCEPT** → this action performs the perception of the pouches obtaining as output their position (x, y) and also the number recognition of any opened drawer
- **SCALE – TIGHTNESS - STRENGTH** → this actions performs the interaction with the instruments for carrying out the corresponding test

 **Movelt**





YOLOv5



# Actions & Fluents

- **PRESENTPOUCH** → returns *true* if at least one pouch has been detected, *false* otherwise
- **SWITCHDRAWER** → returns *true* if at least one drawer is available to be open, *false* otherwise
- **SENSEPOUCH** → returns *true* if the pouch has been grasped, *false* otherwise
- **SENSEOPENBOX** → returns *true* if all the ArUco markers has been detected and the drawer number has been perceived, *false* otherwise
- **SENSEARRUCO** → returns *true* if the number of ArUco markers detected is equal to 0 or 4, *false* otherwise
- **SENSE{SCALE – TIGHTNESS – STRENGTH}** → returns *true* if the value perceived during the test exceeds a threshold, *false* otherwise



# Robot movements overview

BOX open  
GOTO vision2



GOTO vision2  
PERCEPT pouch



BOX open  
GOTO vision2



GOTO pouch  
MOVEGRIPPER  
close



GOTO bin  
MOVEGRPPER  
open




GOTO scale  
GOTO mark10



GOTO mark10\_out  
GOTO mark10\_in



GOTO mark10  
GOTO mark10\_in



# Plan Implementation

**Petri Net Plans (PNP)** → formalism for high level description of plans

**Plan Execution Interface (PLEXI)** → frameworks for implementation of actions/predicates


```
< presentpouch? say_PresentPouch :  
  (not presentpouch)? say_NoPresentPouch;  
  movegripper_o; box_close; goto_vision2;  
  GOTO_LABEL_DRAWER  
>;  
movegripper_20;  
goto_pouch_pose;  
movegripper_155;  
wait_1;  
goto_pickup;  
goto_vision2;  
percept_pouch;  
< sensepouch? say_SensedPouch :  
  (not sensepouch)? say_HelpRequestMissPouch;  
  finalpglabrecovery >;  
scale_reset;  
goto_scale_up;
```

Very  
complex!!!

**Short** sample of the **NOMINAL** plan... from which a **recovery procedure** resume its execution if any action fails

This plan implementation will allow future theses to use a plan generation system that will use the same actions/fluents implemented

# Recovery procedure flow






# Enjoy the video 😊

ALPHA - APPLICATION PG LAB


Execution of PGLab | June 30, 2022 12:37:13

**PAUSE** **STOP** **EXIT**


**GRIPPER CAMERA**



**LAB CAMERA**



**POUCH DETECTED**



**ROBOT PLAN LOG**

POUCH ID: 66 DRAWER SEL: 4 5 6 LINE: [6, 3]

TOTAL\_POUCH: 5  
WORKED\_POUCH: 3  
MISSING\_POUCH: 2  
CURRENT\_DRAWER: 6

goto vision2  
percept pouch


say HelpRequestMissPouch

**DATA ACQUIRED FROM DEVICES**

**Restarted the robot** **SS**

Value: 23.94 (g) Value: 25.21 (mm)

# Thanks for watching!!!



**SAPIENZA**  
UNIVERSITÀ DI ROMA